

# FINGERTIP PULSE OXIMETER

## HF-84A1

# HIFMED®

### Features

- Product adopts double color OLED display , can show the six different display mode.
- Low-power consumption, continuously work for more than six hours with two AAA batteries
- Low voltage indicator
- In the absence of signals, the product will be in after 8 seconds to enter standby state.
- Small in volume, light in weight, and convenient to carry

### SpO2

- Measurement range: 70%~100%
- Accuracy: 80%~100%:±2% (Including 80%);  
70%~79%: ±3% ;Below 70% no requirement;
- Resolution: 1%

### Pulse Rate

- PR: Measurement range: 30BPM~254BPM
- Accuracy: ≤100BPM ,±1BPM;>100BPM ,±2BPM

### Environment&Packaging Information

- Power Requirements:2pcs AAA-size alkaline batteries
- Power: <30mA
- Operation environment:Temperature: 5 C ~ 40 C
- Humidity: 15% ~ 80%,no condensing
- Atmospheric pressure: 70kPa ~ 106kPa
- Transport, storage environment:Temperature: -10 C ~ 40 C
- Humidity: 10%~ 95%,no condensing
- Atmospheric pressure: 70kPa ~ 106kPa,
- Dimension : about 58mmx36mmx33mm
- Weight: 33.7g( without battery)

### Package including

- The host 1 set
- Lanyard 1 pcs
- One user manual 1 pcs
- AAA battery (optional) 2pcs



## FINGER TIP PULSE OXIMETER

SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD.

Add:Building2, NO.885,Qiuxing Road,Nicheng Town, Pudong New Area,Shanghai,China.

Sales department

Tel:+8621-58082099, 58084992

Fax:+8621-58085105

E-mail: lisa@sh-huifeng.com , lisa@huifeng-co.com

Zip:201306

www.sh-huifeng.com

## HF-84A1